

IMPORTANT SAFETY AND PRODUCT INFORMATION

To set up and maintain your VO2 Master Electronic Device (the “VO2 Device”), go to: www.vo2master.com/support

Your VO2 Device comes with the following information booklets:

1. Quick Start Guide,
2. Limited Warranty Information, and
3. This booklet containing Important Safety and Production Information.

CAUTION: For your safety, always read and follow all operating and cleaning instructions carefully before using your VO2 Device.

CAUTION: Consult your doctor before use if you have any pre-existing medical or health conditions that might be affected by your use of this VO2 Device.


WARNING

1. Use this device only for its intended use as described in the Quick Start Guide. Failure to do so may result in serious injury.
2. Do not expose your VO2 Device to extremely high (above 45 degrees Celsius) or low (below 0 degrees Celsius) temperatures.
3. Do not use your VO2 Device in a sauna, steam room, or shower.
4. Do not use your VO2 Device outdoors.
5. Do not place in, or submerge into, water or other liquids.
6. Do not insert any object into, or otherwise block, any opening on the VO2 Device. Keep all openings free of hair, dirt, or other debris.
7. Do not use your VO2 Device if it has been dropped or damaged, submerged into water, or if you see any cracks on the mouthpiece. If such an event occurs, return the device to an authorized service center for examination and repair.
8. Do not use abrasive cleaners to clean your VO2 Device.
9. Remove your VO2 Device if it feels warm or hot.
10. Substances in this product and its battery may harm the environment or cause injury if handled and disposed of improperly.
11. Do not use or operate the VO2 Device while driving or in other situations where distractions could be hazardous. Always be aware of your surroundings when exercising.

12. This product is not a toy. Do not allow children or pets to play with your VO2 Device. The product may contain small components that can be a choking hazard.
13. The VO2 Device may pose risks to users with certain health conditions. Consult your doctor prior to use if you:
 - a. have a heart, respiratory, or other cardiovascular condition;
 - b. are taking any medication;
 - c. suffer from claustrophobia;
 - d. have epilepsy or are sensitive to blinking lights;
 - e. have reduced circulation or bruise easily; or
 - f. have neck, upper back, or other musculoskeletal disorders.
14. Consult your doctor before beginning or modifying any exercise program.
15. Risks of using the VO2 Device may include, but are not limited to:
 - a. inhaled or contacted contaminants (particulates or microbes) causing trauma or illness;
 - b. irritation, abrasion, or ulceration caused by the friction of the facial mask;
 - c. biological effect on mucous, lungs, or skin;
 - d. excessive re-breathing of exhaled air from increased airway dead space or breathing resistance causing an imbalance of blood gases, strain, or atelectasis;
 - e. flow diverted or leaked from or into the oxygen sensor gas circuit causing erroneous test results.
16. If you notice any skin irritation, difficulty breathing, or other health or medical irregularity while or after wearing your VO2 Device, remove your VO2 Device immediately, discontinue use, and consult your doctor.
17. **This product is not a medical device and is not intended to diagnose, treat, cure, or prevent any disease. With regard to accuracy, the VO2 Device has leading hardware and algorithms to track volume of oxygen consumed by the body. We are constantly improving our products to calculate measurements as accurately as possible. The accuracy of VO2 Devices is not intended to match medical devices or scientific measurement devices, but is intended to give you the best information available in a wearable VO2 analyzer.**

Regulatory Information

Federal Communications Commission (FCC) Statement [USA]



Contains FCC ID: 2AA9B04

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- 1) Reorient or relocate the receiving antenna.
- 2) Increase the separation between the equipment and receiver.
- 3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- 4) Consult the dealer or an experienced radio/TV technician for help.

Caution: changes or modifications to the VO2 Device could void the user's authority to operate the equipment.

Industry Canada (IC) Information

Contains IC ID: 12208A-04

This device complies with Industry Canada's license-exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

Declaration of Conformity with Regard to the EU Directive 1999/5/EC

This device complies with EU Directive 1999/5/EC.

For other regional-specific compliance information, please contact us by visiting: www.vo2master.com/support